

FUME HOODS



- **FUME HOODS**

- **Style**

- Bench Hood
- Distillation Hood
- Walk In Type Hood

- **Exhaust Style**

- Constant Volume By-pass
- Restricted By-pass Variable Air Volume
- Auxiliary Air Supply


- **Sash Configuration**

- Standard Vertical Sliding Unframed
- Standard Vertical Sliding Framed
- Horizontal Sliding
- Combination of Vertical and Horizontal Sliding



SPECIAL PURPOSE HOODS

- Perchloric Acid Hoods
- Radioisotope Hoods
- Canopy Hoods
- Spot Extractors



SPOT EXTRACTOR

